

Rabbit Anti-Human Cd317/Bst2 Monoclonal Antibody For Fc - Q10589

Item Number: CM0010677



Introduction

-conjugated rabbit monoclonal antibody against human CD317/BST2 (tetherin), a viral budding restriction factor and multiple myeloma marker. For flow cytometry-based immunophenotyping and antiviral research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD317/BST2 mAb
Remarks/Alias	CD317; HM1.24; TETHERIN
Species	Human
GeneID (Human)	684
GeneID	684
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 49-160 of human CD317/BST2 (NP_004326.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q10589
Protein Weight	18kDa/19kDa
Shipping	Ice bag