

Rabbit Anti-Human Cd31/Pecam1 Mab For Fc - P16284

Item Number: CM0008784



Introduction

-conjugated rabbit monoclonal antibody targeting human CD31/PECAM1 for flow cytometry (FC) applications. Validated on human samples, recognizes an 83 kDa protein. Ideal for leukocyte and endothelial cell studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD31/PECAM1 mAb
Remarks/Alias	CD31; PECA1; GPIIA'; PECAM-1; endoCAM; CD31/EndoCAM
Species	Human
GeneID (Human)	5175
GeneID	5175
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 28-601 of human CD31 (NP_000433.4).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P16284
Protein Weight	83kDa
Shipping	Ice bag