

Rabbit Anti-Human Cd326/Epcam Monoclonal Antibody For Fc - P16422

Item Number: CM0005021



Introduction

-conjugated rabbit monoclonal antibody against human CD326/EPCAM, validated for flow cytometry. Detects epithelial cell adhesion molecule for cancer and stem cell research applications.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD326/EPCAM mAb
Remarks/Alias	ESA; KSA; M4S1; MK-1; DIAR5; EGP-2; EGP40; KS1/4; MIC18; TROP1; BerEp4; EGP314; HNPCC8; LYNCH8; MOC-31; Ber-Ep4; TACSTD1
Species	Human
GeneID (Human)	4072
GeneID	4072
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-265 of human CD326/EPCAM (NP_002345.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P16422
Protein Weight	35kDa
Shipping	Ice bag