

Rabbit Anti-Human Cd33 Monoclonal Antibody For Flow Cytometry - P20138

Item Number: CM0004807



Introduction

-conjugated rabbit monoclonal antibody against human CD33 (Siglec-3, gp67). Recognizes the extracellular domain. Validated for flow cytometry. Ideal for immunophenotyping of myeloid cells and microglia. Research use only.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD33 mAb
Remarks/Alias	CD33; SIGLEC-3; SIGLEC3; p67; CD33 molecule
Species	Human
Gene ID (Human)	945
Gene ID	945
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 18-259 of human CD33(NP_001763.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P20138
Protein Weight	25kDa/33kDa/39kDa
Shipping	Ice bag