

Rabbit Anti-Human Cd335/Nkp46 Monoclonal Antibody For Fc - O76036

Item Number: CM0012537



Introduction

Recombinant rabbit monoclonal antibody targeting human CD335/NKp46, validated for flow cytometry. Recognizes natural cytotoxicity triggering receptor 1 on NK cells. Research use only.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD335/NKp46 mAb
Remarks/Alias	LY94; CD335; NKP46; NK-p46
Species	Human
Gene ID (Human)	9437
Gene ID	9437
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 22-254 of human CD335/NKp46 (NP_004820.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	O76036
Protein Weight	21kDa/22kDa/23kDa/32kDa/34kDa
Shipping	Ice bag