

Rabbit Anti-Human Cd335/Nkp46 Monoclonal Antibody For If/Icc And Fc - O76036

Item Number: CM0009028



Introduction

-conjugated rabbit monoclonal antibody against human CD335/NKp46 (NCR1). Suitable for IF/ICC and flow cytometry detection of NK cells. Validated in human samples.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD335/NKp46 mAb
Remarks/Alias	NCR1; CD335; LY94; NK-p46; NKP46; natural cytotoxicity triggering receptor 1
Species	Human
GeneID (Human)	9437
GeneID	9437
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-254 of human CD335/NKp46 (NP_004820.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC
Cross Reactivity	Human
SWISS	O76036
Protein Weight	21kDa/22kDa/23kDa/32kDa/34kDa
Shipping	Ice bag