

Rabbit Anti-Human Cd335/Nkp46 Monoclonal Antibody For If/Icc, Fc - O76036

Item Number: CM0010427



Introduction

-conjugated rabbit monoclonal antibody against human CD335/NKp46 (NCR1), validated for immunofluorescence and flow cytometry. Ideal for natural killer cell research and immunophenotyping.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD335/NKp46 mAb
Remarks/Alias	NCR1, CD335, LY94, NK-p46, NKP46, natural cytotoxicity triggering receptor 1
Species	Human
Gene ID (Human)	9437
Immunogen	Recombinant protein corresponding to amino acids 22-254 of human CD335/NKp46 (NP_004820.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC
Cross Reactivity	Human
SWISS	O76036
Protein Weight	21kDa/22kDa/23kDa/32kDa/34kDa (multiple isoforms)
Shipping	Ice bag