

Rabbit Anti-Human Cd336/Nk-P44 Mab For Fc - O95944

Item Number: CM0010171



Introduction

High-quality -conjugated rabbit monoclonal antibody against human CD336/NK-p44 (NCR2). Validated for flow cytometry (FC) of activated natural killer cells; 31kDa. Ideal for tumor immunology and NK cell research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD336/NK-p44 mAb
Remarks/Alias	LY95; CD336; NKP44; NK-p44; dj149M18.1
Species	Human
GeneID (Human)	9436
GeneID	9436
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-190 of human CD336/NK-p44 (NP_004819.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	O95944
Protein Weight	31kDa
Shipping	Ice bag