

Rabbit Anti-Human Cd337/Nkp30 Monoclonal Antibody For Flow Cytometry - O14931

Item Number: CM0004132



Introduction

High-quality rabbit monoclonal antibody targeting human CD337/NKp30 (NCR3) for flow cytometry. Ideal for studying NK cell cytotoxicity, tumor immunology, and dendritic cell maturation.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD337/NKp30 mAb
Remarks/Alias	1C7; MALS; CD337; LY117; NKp30
Species	Human
GeneID (Human)	259197
GeneID	259197
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 19-135 of human CD337/NKp30 (NP_667341.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (Flow Cytometry)
Cross Reactivity	Human
SWISS	O14931
Protein Weight	16kDa/17kDa/18kDa/19kDa/20kDa/21kDa
Shipping	Ice bag