

Rabbit Anti-Human Cd337/Nkp30 Monoclonal Antibody For If/Icc And Fc (Intra) - O14931

Item Number: CM0006067



Introduction

-conjugated rabbit anti-human CD337/NKp30 monoclonal antibody validated for IF/ICC and intracellular flow cytometry. Targets natural cytotoxicity triggering receptor 3 (NCR3) expressed on NK cells; human reactive, O14931.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD337/NKp30 mAb
Remarks/Alias	NCR3; 1C7; CD337; LY117; MALS; NKp30; natural cytotoxicity triggering receptor 3
Species	Human
GeneID (Human)	259197
GeneID	259197
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-135 of human CD337/NKp30 (NP_667341.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC (intra)
Cross Reactivity	Human
SWISS	O14931
Protein Weight	16kDa/17kDa/18kDa/19kDa/20kDa/21kDa
Shipping	Ice bag