

Rabbit Anti-Human Cd34 Monoclonal Antibody For Fc - P28906

Item Number: CM0009091

Introduction

-conjugated rabbit monoclonal antibody against human CD34 (P28906) for flow cytometry. Detects hematopoietic progenitor cells and endothelial cells. Suitable for FC applications in hematopoietic research. Research use only.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD34 mAb
Remarks/Alias	CD34;CD34 molecule;GIG3
Species	Human
Gene ID (Human)	947
Gene ID	947
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 170-290 of human CD34 (NP_001020280.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P28906
Protein Weight	35kDa/41kDa
Shipping	Ice bag