

Rabbit Anti-Human Cd34 Monoclonal Antibody For If/Icc And Fc - P28906

Item Number: CM0004865

Introduction

-labeled rabbit monoclonal antibody targeting human CD34, validated for immunofluorescence (IF/ICC) and flow cytometry (FC). Ideal for hematopoietic stem cell and endothelial cell studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD34 mAb
Remarks/Alias	CD34;CD34 molecule;GIG3;MORT1
Species	Human
Gene ID (Human)	947
Gene ID	947
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 170-290 of human CD34 (NP_001020280.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF/ICC, FC
Cross Reactivity	Human
SWISS	P28906
Protein Weight	35kDa/41kDa
Shipping	Ice bag