

Rabbit Anti-Human Cd35 Mab For Fc - P17927

Item Number: CM0005027

Introduction

Rabbit anti-Human CD35 monoclonal antibody for flow cytometry detection of human complement receptor type 1 (CR1) on erythrocytes and immune cells. Research-use only.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD35 mAb
Remarks/Alias	CR1; C3BR; C4BR; CD35; KN; complement receptor type 1
Species	Human
GeneID (human)	1378
GeneID	1378
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 901-1095 of human CD35 (NP_000564.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P17927
Protein Weight	223kDa
Shipping	Ice bag