

## Rabbit Anti-Human Cd35 Monoclonal Antibody For Fc - P17927

Item Number: CM0007933



### Introduction

-conjugated rabbit monoclonal antibody targeting human CD35 (CR1). Validated for flow cytometry. Ideal for studying complement-mediated immune adherence, opsonization, viral entry, and regulatory T cell differentiation.

[Learn More](#)

| Parameter        | Value                                                                                              |
|------------------|----------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Human CD35 mAb                                                                         |
| Remarks/Alias    | CR1; C3BR; C4BR; CD35; KN; complement receptor type 1                                              |
| Species          | Human                                                                                              |
| Gene ID (Human)  | 1378                                                                                               |
| Gene ID          | 1378                                                                                               |
| Immunogen        | Recombinant protein (recombinant fusion protein containing aa 901-1095 of human CD35, NP_000564.2) |
| Source           | Rabbit                                                                                             |
| Category         | Monoclonal Antibodies                                                                              |
| Application      | FC                                                                                                 |
| Cross Reactivity | Human                                                                                              |
| SWISS            | P17927                                                                                             |
| Protein Weight   | 223kDa                                                                                             |
| Shipping         | Ice bag                                                                                            |