

Rabbit Anti-Human Cd35 Monoclonal Antibody For Flow Cytometry - P17927

Item Number: CM0008365



Introduction

-conjugated rabbit monoclonal antibody targeting human CD35 (CR1). Validated for flow cytometry, it detects complement receptor type 1 on erythrocytes and immune cells. Ideal for studying immune adherence and complement-mediated clearance.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD35 mAb
Remarks/Alias	CR1; C3BR; C4BR; CD35; KN; complement receptor type 1
Species	Human
Gene ID (Human)	1378
Gene ID	1378
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 901-1095 of human CD35 (NP_000564.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P17927
Protein Weight	223kDa
Shipping	Ice bag