


Rabbit Anti-Human Cd352/Slamf6 Monoclonal Antibody For Fc - Q96Du3

Item Number: CM0012043



Introduction

-conjugated rabbit monoclonal antibody targeting human CD352/SLAMF6 for flow cytometry. Validated for human samples; SLAMF6 is a self-ligand receptor modulating NK and T-cell activation.

[Learn More](#)

Parameter	Specification
Product Name	Rabbit anti-Human CD352/SLAMF6 mAb
Remarks/Alias	SLAMF6; CD352; KALI; KALib; Ly108; NTB-A; NTBA; SF2000; SLAM family member 6
Species	Human
GeneID (Human)	114836
GeneID	114836
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-215 of human CD352/SLAMF6 (NP_001171643.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q96DU3
Protein Weight	24kDa/37kDa
Shipping	Ice bag