

Rabbit Anti-Human Cd353/Slamf8 Monoclonal Antibody For Flow Cytometry - Q9P0V8

Item Number: CM0009379

Introduction

-conjugated rabbit monoclonal antibody targeting human CD353 (SLAMF8). Validated for flow cytometry, suitable for B-cell and lymphoid tissue research. Reliable IVD raw material from CamelBio.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD353/SLAMF8 mAb
Remarks/Alias	BLAME; CD353; SBB142
Species	Human
GeneID (Human)	56833
GeneID	56833
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-123 of human CD353/SLAMF8 (NP_064510.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9P0V8
Protein Weight	32kDa
Shipping	Ice bag