

Rabbit Anti-Human Cd369/Dectin-1 Monoclonal Antibody For Flow Cytometry - Q9Bxn2

Item Number: CM0008914

Introduction

-conjugated rabbit monoclonal antibody against human CD369/Dectin-1 (CLEC7A) for flow cytometry. Targets the extracellular domain; ideal for detecting Dectin-1 on leukocytes and dendritic cells in innate immunity and fungal infection research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD369/Dectin-1 mAb
Remarks/Alias	CLEC7A; BGR; CANDF4; CD369; CLECSF12; DECTIN1; SCARE2; C-type lectin domain containing 7A
Species	Human
GeneID (human)	64581
GeneID	64581
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 66-247 of human CD369/Dectin-1 (NP_922938.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q9BXN2
Protein Weight	5kDa/8kDa/13-27kDa
Shipping	Ice bag