


## Rabbit Anti-Human Cd370/Clec9A Mab - Q6Uxn8

Item Number: CM0007863



### Introduction

Highly specific conjugated rabbit monoclonal antibody targeting human CLEC9A/CD370, validated for flow cytometry. Ideal for dendritic cell and dead-cell antigen cross-presentation research.

[Learn More](#)

| Parameter        | Value                                                                                                                     |
|------------------|---------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Human CD370/CLEC9A mAb                                                                                        |
| Remarks/Alias    | CD370; DNGR1; DNGR-1; UNQ9341                                                                                             |
| Species          | Human                                                                                                                     |
| Gene ID          | 283420                                                                                                                    |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 110-241 of human CD370/CLEC9A (NP_997228.1) |
| Source           | Rabbit                                                                                                                    |
| Category         | Monoclonal Antibodies                                                                                                     |
| Application      | FC                                                                                                                        |
| Cross Reactivity | Human                                                                                                                     |
| SWISS            | Q6UXN8                                                                                                                    |
| Protein Weight   | 27kDa                                                                                                                     |
| Shipping         | Ice bag                                                                                                                   |