


Rabbit Anti-Human Cd370/Clec9A Monoclonal Antibody For Flow Cytometry - Q6Uxn8

Item Number: CM0011780



Introduction

-conjugated rabbit monoclonal antibody against human CD370/CLEC9A (DNGR-1), validated for flow cytometry. Recognizes a dendritic cell receptor for dead cell-associated antigens. Ideal for immunophenotyping and biological research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD370/CLEC9A mAb
Remarks/Alias	CD370; DNGR1; DNGR-1; UNQ9341
Species	Human
GeneID (Human)	283420
GeneID	283420
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 110-241 of human CD370/CLEC9A (NP_997228.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q6UXN8
Protein Weight	27kDa
Shipping	Ice bag