


Rabbit Anti-Human Cd370/Clec9A Monoclonal Antibody For Flow Cytometry - Q6Uxn8

Item Number: CM0006042



Introduction

-conjugated rabbit monoclonal antibody targeting human CD370/CLEC9A (Q6UXN8). Validated for flow cytometry, enabling identification of BDCA3+ dendritic cells and studies of dead cell antigen cross-presentation.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD370/CLEC9A mAb
Remarks/Alias	CD370; DNGR1; DNGR-1; UNQ9341
Species	Human
GeneID	283420
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 110-241 of human CD370/CLEC9A (NP_997228.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q6UXN8
Protein Weight	27kDa
Shipping	Ice bag