


Rabbit Anti-Human Cd38 Mab - P28907

Item Number: CM0009353



Introduction

Rabbit anti-Human CD38 monoclonal antibody for flow cytometry. Targets CD38 (ADP-ribosyl cyclase 1), involved in calcium signaling and cell adhesion. Suitable for human samples.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD38 mAb
Remarks/Alias	CD38; ADPRC 1; ADPRC1; CD38 molecule
Species	Human
Gene ID (Human)	952
Gene ID	952
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 43-300 of human CD38 (P28907)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P28907
Protein Weight	13kDa/34kDa
Shipping	Ice bag