


Rabbit Anti-Human Cd38 Monoclonal Antibody For Flow Cytometry - P28907

Item Number: CM0005572



Introduction

-conjugated rabbit monoclonal antibody targeting human CD38 for flow cytometry. CD38 is a multifunctional transmembrane glycoprotein with enzymatic and receptor functions, involved in calcium signaling and lymphocyte activation. Ideal for immunophenotyping and research on immune cells, diabetes, and cancer.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD38 mAb
Remarks/Alias	CD38; ADPRC 1; ADPRC1; CD38 molecule
Species	Human
GeneID (Human)	952
GeneID	952
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 43-300 of human CD38 (P28907).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P28907
Protein Weight	13kDa/34kDa
Shipping	Ice bag