

Rabbit Anti-Human Cd39/Entpd1 Mab For Fc - P49961

Item Number: CM0008818

Introduction

-conjugated rabbit monoclonal antibody targeting human CD39/ENTPD1 for flow cytometry. Specific for human CD39, validated for FC applications. Ideal for immunological studies of ectonucleotidase activity and lymphocyte activation.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD39/ENTPD1 mAb
Remarks/Alias	CD39; SPG64; ATPDase; NTPDase-1
Species	Human
GeneID (Human)	953
GeneID	953
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 255-442 of human CD39/ENTPD1 (NP_001767.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P49961
Protein Weight	58kDa
Shipping	Ice bag