

Rabbit Anti-Human Cd39/Entpd1 Mab For Fc - P49961

Item Number: CM0004899

Introduction

-conjugated rabbit monoclonal antibody targeting human CD39/ENTPD1 for flow cytometry. Validated in human samples for ectonucleotidase and immune activation research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD39/ENTPD1 mAb
Remarks/Alias	CD39; SPG64; ATPDase; NTPDase-1
Species	Human
GeneID	953 (Human)
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 255-442 of human CD39/ENTPD1 (NP_001767.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (Flow Cytometry)
Cross Reactivity	Human
SWISS	P49961
Protein Weight	58kDa
Shipping	Ice bag