


Rabbit Anti-Human Cd40 Mab - P25942

Item Number: CM0008514



Introduction

Rabbit anti-Human CD40 monoclonal antibody for flow cytometry (FC). Detects human CD40 (p50, Bp50, CDW40, TNFRSF5), a key immune receptor. Ideal for immunology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD40 mAb
Remarks/Alias	p50; Bp50; CDW40; TNFRSF5
Species	Human
GeneID (human)	958
GeneID	958
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-193 of human CD40 (P25942).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P25942
Protein Weight	31kDa
Shipping	Ice bag