


Rabbit Anti-Human Cd40 Monoclonal Antibody For Fc - P25942

Item Number: CM0011134



Introduction

-conjugated rabbit monoclonal antibody against human CD40, validated for flow cytometry. Targets B-cell surface antigen CD40 (Bp50, CDw40, TNFRSF5) with human cross-reactivity. Ideal for direct detection in immunology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD40 mAb
Remarks/Alias	p50; Bp50; CDW40; TNFRSF5
Species	Human
Gene ID (Human)	958
Gene ID	958
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-193 of human CD40 (NP_001241.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P25942
Protein Weight	31kDa
Shipping	Ice bag