


Rabbit Anti-Human Cd40 Monoclonal Antibody For Flow Cytometry - P25942

Item Number: CM0012093



Introduction

-conjugated rabbit monoclonal antibody targeting human CD40 (TNFRSF5), validated for flow cytometry. Ideal for studying B-cell and macrophage activation, immune response, and CD40-CD40L signaling.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD40 mAb
Remarks/Alias	p50; Bp50; CDW40; TNFRSF5
Species	Human
GeneID (Human)	958
GeneID	958
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-193 of human CD40 (P25942)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P25942
Protein Weight	31kDa
Shipping	Ice bag