

Rabbit Anti-Human Cd42B Monoclonal Antibody For Flow Cytometry - P07359

Item Number: CM0011433



Introduction

-conjugated rabbit anti-human CD42b monoclonal antibody for flow cytometry (FC). Targets platelet glycoprotein Ib alpha chain (GP1BA, CD42b, P07359). Suitable for hemostasis, thrombosis, and platelet biology research. High-quality IVD raw material.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD42b mAb
Remarks/Alias	BSS; GP1B; VWDP; CD42B; GPIbA; BDPLT1; BDPLT3; DBPLT3; GPIbalpha; CD42b-alpha
Species	Human
GeneID (Human)	2811
GeneID	2811
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 15-531 of human CD42b (NP_000164.5)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P07359
Protein Weight	72kDa
Shipping	Ice bag