

Rabbit Anti-Human Cd45Ra Monoclonal Antibody () For Flow Cytometry - P08575-8

Item Number: CM0012168

Introduction

-conjugated rabbit monoclonal antibody against human CD45RA (isoform 6) for flow cytometry, 138kDa. Ideal for human leukocyte subset identification and immunophenotyping.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD45RA mAb
Remarks/Alias	LCA; LY5; B220; CD45; L-CA; T200; CD45R; GP180; IMD105
Species	Human
GeneID (human)	5788
GeneID	5788
Immunogen	Recombinant protein (Recombinant fusion protein containing a sequence corresponding to amino acids 26-482 of human CD45RA, P08575-8)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P08575-8
Protein Weight	138kDa
Shipping	Ice bag