

Rabbit Anti-Human Cd45Ra Monoclonal Antibody For Flow Cytometry - P08575-8

Item Number: CM0012206

Introduction

-conjugated rabbit anti-human CD45RA monoclonal antibody for flow cytometry. Detects the 138 kDa CD45 isoform 6 (L-CA, T200) on human leukocytes. Reliable tool for immunophenotyping and immune cell analysis.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD45RA mAb
Remarks/Alias	LCA; LY5; B220; CD45; L-CA; T200; CD45R; GP180; IMD105
Species	Human
GenelD (Human)	5788
GenelD	5788
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-482 of human CD45RA (P08575-8)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P08575-8
Protein Weight	138kDa
Shipping	Ice bag