

Rabbit Anti-Human Cd45Ra Monoclonal Antibody For Flow Cytometry - P08575-8

Item Number: CM0008797

Introduction

High-quality -conjugated rabbit monoclonal antibody targeting human CD45RA (PTPRC isoform 6). Validated for flow cytometry, enabling precise leukocyte immunophenotyping in clinical and research settings.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD45RA mAb
Remarks/Alias	LCA; LY5; B220; CD45; L-CA; T200; CD45R; GP180; IMD105
Species	Human
Gene ID (Human)	5788
Gene ID	5788
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-482 of human CD45RA (P08575-8).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P08575-8
Protein Weight	138kDa
Shipping	Ice bag