

Rabbit Anti-Human Cd45Ra Monoclonal Antibody For Flow Cytometry - P08575-8

Item Number: CM0009781

Introduction

-conjugated rabbit monoclonal antibody targeting human CD45RA (isoform 6, P08575-8) for flow cytometry. Cross-reacts with human; ideal for leukocyte immunophenotyping studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD45RA mAb
Remarks/Alias	LCA; LY5; B220; CD45; L-CA; T200; CD45R; GP180; IMD105
Species	Human
Gene ID (human)	5788
Gene ID	5788
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-482 of human CD45RA (P08575-8)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P08575-8
Protein Weight	138kDa
Shipping	Ice bag