

## Rabbit Anti-Human Cd45Ra Monoclonal Antibody For Flow Cytometry - P08575-8

Item Number: CM0005057

### Introduction

-conjugated rabbit monoclonal antibody targeting human CD45RA (isoform 6 of PTPRC, P08575-8) for flow cytometry. Ideal for immunophenotyping of hematopoietic cells.

[Learn More](#)

| Parameter        | Value                                                                                                                                |
|------------------|--------------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Human CD45RA mAb                                                                                                         |
| Remarks/Alias    | LCA; LY5; B220; CD45; L-CA; T200; CD45R; GP180; IMD105                                                                               |
| Species          | Human                                                                                                                                |
| GeneID (Human)   | 5788                                                                                                                                 |
| GeneID           | 5788                                                                                                                                 |
| Immunogen        | Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 26-482 of human CD45RA (P08575-8). |
| Source           | Rabbit                                                                                                                               |
| Category         | Monoclonal Antibodies                                                                                                                |
| Application      | FC                                                                                                                                   |
| Cross Reactivity | Human                                                                                                                                |
| SWISS            | P08575-8                                                                                                                             |
| Protein Weight   | 138kDa                                                                                                                               |
| Shipping         | Ice bag                                                                                                                              |