

Rabbit Anti-Human Cd45Rb Monoclonal Antibody For Fc - P08575-9

Item Number: CM0009555

Introduction

-conjugated rabbit monoclonal antibody targeting human CD45RB (isoform 7 of PTPRC), validated for flow cytometry. Ideal for leukocyte common antigen research and immune cell profiling.

[Learn More](#)

Parameter	Specification
Product Name	Rabbit anti-Human CD45RB mAb
Remarks/Alias	LCA; LY5; B220; CD45; L-CA; T200; CD45R; GP180; IMD105
Species	Human
GeneID (Human)	5788
GeneID	5788
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P08575-9
Shipping	Ice bag