

Rabbit Anti-Human Cd45Rb Monoclonal Antibody (Pe) For Flow Cytometry - P08575-9

Item Number: CM0010888

Introduction

PE-conjugated rabbit monoclonal antibody targeting human CD45RB (PTPRC isoform 7, SWISS: P08575-9), for flow cytometry. Recombinant immunogen (aa 26-463), cross-reacts with human. ~136 kDa. Suitable for leukocyte research.

[Learn More](#)

Parameter	Value
Product Name	PE Rabbit anti-Human CD45RB mAb
Remarks/Alias	LCA; LY5; B220; CD45; L-CA; T200; CD45R; GP180; IMD105
Species	Human
GeneID	5788
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 26-463 of human CD45RB (XP_006711537.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P08575-9
Protein Weight	136kDa
Shipping	Ice bag