


Rabbit Anti-Human Cd46 Mab For Flow Cytometry - P15529

Item Number: CM0009422



Introduction

-conjugated rabbit monoclonal antibody targeting human CD46 (MCP). Validated for flow cytometry, this antibody recognizes human CD46, a key complement regulator and pathogen receptor. Ideal for immunology and infectious disease research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD46 mAb
Remarks/Alias	MCP; TLX; AHUS2; MIC10; TRA2.10
Species	Human
GeneID (human)	4179
GeneID	4179
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-328 of human CD46 (NP_002380.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P15529
Protein Weight	37-44kDa
Shipping	Ice bag