


Rabbit Anti-Human Cd46 Monoclonal Antibody For Fc - P15529

Item Number: CM0011838



Introduction

-conjugated rabbit monoclonal antibody targeting human CD46 (MCP). Validated for flow cytometry; recognizes a 37-44 kDa membrane protein. Suitable for complement research, viral entry studies, and immune regulation assays.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD46 mAb
Remarks/Alias	MCP; TLX; AHUS2; MIC10; TRA2.10
Species	Human
Gene ID (Human)	4179
Gene ID	4179
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-328 of human CD46 (NP_002380.3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P15529
Protein Weight	37-44kDa
Shipping	Ice bag