

Rabbit Anti-Human Cd47 Monoclonal Antibody For Fc - Q08722

Item Number: CM0011656



Introduction

-conjugated rabbit monoclonal antibody targeting human CD47, validated for flow cytometry. Supports immunology and IVD assay development; cross-reacts with human. From CamelBio's raw-material portfolio. (SWISS: Q08722)

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD47 mAb
Remarks/Alias	CD47; IAP; MER6; OA3; CD47 molecule
Species	Human
GeneID	961
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-139 of human CD47 (NP_001768.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q08722
Protein Weight	31kDa/33kDa/35kDa
Shipping	Ice bag