


Rabbit Anti-Human Cd47 Monoclonal Antibody For Flow Cytometry - Q08722

Item Number: CM0008498



Introduction

-conjugated rabbit anti-human CD47 monoclonal antibody for flow cytometry. Specific for human CD47 (SWISS: Q08722). Suitable for investigating CD47's role in immune regulation, cancer biology, and cardiovascular function.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD47 mAb
Remarks/Alias	CD47; IAP; MER6; OA3; CD47 molecule
Species	Human
GeneID (Human)	961
GeneID	961
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-139 of human CD47 (NP_001768.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	Q08722
Protein Weight	31kDa/33kDa/35kDa
Shipping	Ice bag