

Rabbit Anti-Human Cd48 Monoclonal Antibody For Fc - P09326

Item Number: CM0008994



Introduction

This -conjugated rabbit monoclonal antibody recognizes human CD48 (BCM1, BLAST, SLAMF2), a GPI-anchored glycoprotein widely expressed on hematopoietic cells. It is suitable for flow cytometry (FC) applications.

Immunogen: recombinant fusion protein (aa 27-220). Ideal for T-cell signaling and immune cell activation studies.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD48 mAb
Remarks/Alias	BCM1; BLAST; hCD48; mCD48; BLAST1; SLAMF2; MEM-102
Species	Human
GeneID (Human)	962
GeneID	962
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-220 of human CD48 (NP_001769.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P09326
Protein Weight	19kDa/28kDa
Shipping	Ice bag