

Rabbit Anti-Human Cd48 Monoclonal Antibody For Flow Cytometry - P09326

Item Number: CM0005517



Introduction

-conjugated rabbit monoclonal antibody targeting human CD48 (SLAMF2). Validated for flow cytometry (FC), this reagent enables reliable detection of CD48 on hematopoietic cells for immunological research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD48 mAb
Remarks/Alias	BCM1; BLAST; hCD48; mCD48; BLAST1; SLAMF2; MEM-102
Species	Human
GeneID (Human)	962
GeneID	962
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-220 of human CD48 (NP_001769.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P09326
Protein Weight	28kDa
Shipping	Ice bag