


Rabbit Anti-Human Cd5 Mab - P06127

Item Number: CM0012495



Introduction

Rabbit anti-Human CD5 monoclonal antibody for flow cytometry (FC). Specific to human CD5 (SWISS: P06127), a T-cell surface glycoprotein involved in immune regulation. Ideal for immunology research.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD5 mAb
Remarks/Alias	T1; LEU1
Species	Human
GeneID (Human)	921
GeneID	921
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-371 of human CD5 (NP_055022.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P06127
Protein Weight	55kDa
Shipping	Ice bag