

Rabbit Anti-Human Cd50(Icam-3) Monoclonal Antibody For Fc - P32942

Item Number: CM0010990



Introduction

-conjugated rabbit monoclonal antibody against human CD50/ICAM-3, optimized for flow cytometry, ideal for studying leukocyte adhesion and dendritic cell interactions.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD50(ICAM-3) mAb
Remarks/Alias	CD50; CDW50; ICAM-R
Species	Human
Gene ID (Human)	3385
Gene ID	3385
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 30-114 of human CD50(ICAM-3) (NP_002153.2)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P32942
Protein Weight	60kDa
Shipping	Ice bag