

Rabbit Anti-Human Cd51 Monoclonal Antibody For Flow Cytometry - P06756

Item Number: CM0012280



Introduction

-conjugated rabbit monoclonal antibody against human CD51 (integrin alpha-V, P06756). Validated for flow cytometry (intracellular). Ideal for research on cell adhesion, integrin signaling, and viral infection mechanisms.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD51 mAb
Remarks/Alias	CD51; MSK8; VNRA; VTNR
Species	Human
GeneID (Human)	3685
GeneID	3685
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 949-1048 of human CD51 (NP_002201.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC (intra)
Cross Reactivity	Human
SWISS	P06756
Protein Weight	116kDa
Shipping	Ice bag