


## Rabbit Anti-Human Cd52 Monoclonal Antibody For Flow Cytometry - P31358

Item Number: CM0012434



### Introduction

High-quality -conjugated rabbit anti-human CD52 monoclonal antibody for flow cytometry. Specific for human CD52 (SWISS P31358), validated for FC. Ideal for immunophenotyping and research applications.

[Learn More](#)

| Parameter        | Value                                                                                             |
|------------------|---------------------------------------------------------------------------------------------------|
| Product Name     | Rabbit anti-Human CD52 mAb                                                                        |
| Remarks/Alias    | HE5; CDW52; EDDM5                                                                                 |
| Species          | Human                                                                                             |
| Gene ID (Human)  | 1043                                                                                              |
| Gene ID          | 1043                                                                                              |
| Immunogen        | Synthetic peptide corresponding to a sequence within amino acids 1-61 of human CD52 (NP_001794.2) |
| Source           | Rabbit                                                                                            |
| Category         | Monoclonal Antibodies                                                                             |
| Application      | FC                                                                                                |
| Cross Reactivity | Human                                                                                             |
| SWISS            | P31358                                                                                            |
| Protein Weight   | 7kDa                                                                                              |
| Shipping         | Ice bag                                                                                           |