


Rabbit Anti-Human Cd52 Monoclonal Antibody For Flow Cytometry - P31358

Item Number: CM0009301



Introduction

-conjugated rabbit monoclonal antibody targeting human CD52 (CAMPATH-1 antigen, HE5). Validated for flow cytometry, specific to human. Swiss-Prot: P31358, ~7 kDa GPI-anchored glycoprotein.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD52 mAb
Remarks/Alias	HE5; CDW52; EDDM5
Species	Human
GeneID (Human)	1043
GeneID	1043
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 1-61 of human CD52 (NP_001794.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P31358
Protein Weight	7kDa
Shipping	Ice bag