


Rabbit Anti-Human Cd55 Monoclonal Antibody For Fc - P08174

Item Number: CM0008685



Introduction

-conjugated rabbit monoclonal antibody against human CD55 (DAF) validated for flow cytometry. Recognizes the GPI-anchored complement regulatory protein; ideal for studies on complement pathways, viral entry, and immune evasion.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD55 mAb
Remarks/Alias	CR; TC; DAF; CROM; CHAPLE
Species	Human
GeneID (Human)	1604
GeneID	1604
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-353 of human CD55 (NP_000565.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P08174
Protein Weight	41kDa
Shipping	Ice bag