


Rabbit Anti-Human Cd55 Monoclonal Antibody For Fc - P08174

Item Number: CM0004517



Introduction

-conjugated rabbit monoclonal antibody targeting human CD55, validated for flow cytometry. Suitable for studying complement regulation, innate immunity, and viral entry. Reacts with human samples and detects multiple CD55 isoforms.

[Learn More](#)

Parameter	Value
Product Name	Rabbit anti-Human CD55 mAb
Remarks/Alias	CR; TC; DAF; CROM; CHAPLE
Species	Human
Gene ID (Human)	1604
Gene ID	1604
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-353 of human CD55 (NP_000565.1).
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P08174
Protein Weight	39kDa/41kDa/48kDa/56kDa/59kDa
Shipping	Ice bag