


Rabbit Anti-Human Cd55 Monoclonal Antibody For Flow Cytometry - P08174

Item Number: CM0008629



Introduction

High-quality -conjugated rabbit monoclonal antibody targeting human CD55 (DAF), validated for flow cytometry. Ideal for complement cascade and host-pathogen interaction studies. – P08174

[Learn More](#)

Parameter	Details
Product Name	Rabbit anti-Human CD55 mAb
Remarks/Alias	CR; TC; DAF; CROM; CHAPLE; CD55 molecule
Species	Human
GeneID (human)	1604
GeneID	1604
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 35-353 of human CD55 (NP_000565.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	FC
Cross Reactivity	Human
SWISS	P08174
Protein Weight	39kDa/41kDa/48kDa/56kDa/59kDa
Shipping	Ice bag